



Thermo Scientific Orion Star A320, A220 and A120-Series Portable Meters

Outside the laboratory, testing brings enough of its own challenges from getting the sample to taking the reading. Just because you test outside the laboratory doesn't mean your testing accuracy or quality needs to suffer. The new Orion Star A series of portable meters features a variety of advantages to handle testing in difficult environments.



Protective plugs cover RS232 and USB ports to prevent corrosion

Bright, backlit LCD display clearly shows results on those bright, sunny days as well as the dark, rainy ones and in dark areas of the plant

IP-67 waterproof housing

Ergonomic design allows for one-handed operation and features a recessed groove in the back of the meter to fit fingers easily wrapping around the meter

Keypad texture works well with gloved and wet hands

Advanced meter design prevents water ingress through all connector ports.

Do you have access to an AC power source while testing?
Save battery life by running the meter off AC power with the optional adapter

