



## Introducing the latest Thermo Scientific Orion Star Meters

Orion VERSA STAR and Star A Meters

The new Orion Star™ A and VERSA STAR™ meters are built on the quality and expertise you've come to expect from the Orion product line. The expanded offering combines the tried-and-true features relied on by many professionals with new features concentrated on making them easier to use and minimizing the time needed to train others or to refer to instructions.



Orion  
VERSA STAR  
Benchtop  
Series

The Orion VERSA STAR advanced benchtop meter and modules exemplify versatility. Four channels can be configured in any combination. Define the channels as you need them. Some suggestions are: four different parameters for complete multiparameter analysis, two sets of two different parameters for simultaneous dual-parameter testing or even four of the same for really high throughput. The versatility also shows in the large, customizable color display; ability to handle the testing spectrum from simple readings to advanced ISE applications; redesigned electrode holder head and exclusive auto-recognition capability for Orion optical RDO and polarographic DO probes. Orion VERSA STAR meters have all the best features and options to work with any application, even the unexpected.

To complement the new meters, we've kept our specialized, industry-proven meters. The Orion DUAL STAR meters are easy-to-use with two channels for simultaneous measurement. For environments that need non-sparking meters, the Orion Intrinsically-Safe Meters meet that criteria and are US intrinsically-safe rated.

Our 2012 Orion meter offering means no more settling for a meter – get the right fit for your application and testing needs.

Orion Star A320 and A200-series meters offer dependability and are an easy choice if you need a dedicated single or multiparameter meter for accurate and repeatable measurements every day. Built with the testing location in mind, the benchtop meters feature an Orion Star stirrer probe port and include the redesigned electrode holder head. The portable meters are waterproof with an IP67-rated housing and have the option to accept AC power. The large, backlit graphic display shows results with plain language prompts. Exclusive auto-recognition capability for Orion optical RDO and polarographic DO probes lets you choose the best probe for your application, not limiting your probe options. These meters are the ideal choice for any lab, for every sample, every day.


Orion Star  
A320 &  
A200 Series











Orion Star  
A100 Series



For simplicity and affordability without sacrificing accuracy, look no further than the Orion Star A100-series meters. These single parameter meters are available in benchtop and portable housings. The simple keypad layout combined with the dedicated onscreen messages help with calibration and setup. Orion Star A100-series meters provide just the core information you need.



Model (Benchtop Meters)	Orion VERSA STAR Premium	Orion DUAL STAR Advanced	Orion Star A200 Series Advanced	Orion Star A100 Series Basic	
Example Display					
Display Type	Customizable color graphic LCD with step-by-step text prompts	Graphic LCD with step-by-step text prompts	Graphic LCD with step-by-step text prompts	Fixed segment LCD with basic text prompts	
Keypad	Numeric with menu-specific function keys	Numeric with menu-specific function keys	Comprehensive with menu-specific function keys	Simple with highlighted measure key	
Number of Channels	1 to 4	2	1 to 2	1	
Orion Star Stirrer Probe Control	2	1	1	—	
Datalog	2000 with date/time	1000 with date/time	2000 with date/time	50	
USB Output	2	1	1	—	
RS232 Output	1	1	1	—	
Methods	10 per channel	10 per channel	10	—	
Parameter	pH	Versa Star 10 meter with pH module Versa Star 80 meter with pH/LogR module Versa Star 50 meter with pH, conductivity modules	Dual Star meter with two pH/ISE channels	Star A211 pH meter Star A215 pH/ conductivity meter Star A216 pH/RDO/DO meter	Star A111 pH meter
	ISE	Versa Star 40A meter with pH/ISE module Versa Star 40B meter with two pH/ISE modules	Dual Star meter with two pH/ISE channels	Star A214 pH/ISE meter	—
	Conductivity/TDS/ Resistivity/Salinity	Versa Star 20 meter with conductivity module Versa Star 50 meter with pH, conductivity modules	—	Star A212 conductivity meter Star A215 pH/ conductivity meter	—
	Conductivity, TDS	—	—	—	Star A112 conductivity meter
	RDO Optical & Polarographic Dissolved Oxygen	Versa Star 30 meter with RDO/DO module	—	Star A213 RDO/DO meter Star A216 pH/RDO/DO meter	—
	Polarographic Dissolved Oxygen	—	—	—	Star A113 DO meter
	pH/ISE/Conductivity/ RDO/DO	Versa Star 90 meter with pH/ISE, conductivity, RDO/DO modules Versa Star 90 meter with two pH/ISE, conductivity, RDO/DO modules	—	—	—
Selectable Meter Components	Versa Star meter customizable with up to four pH modules • pH • pH/ISE • pH/LogR • Conductivity • RDO/DO	—	—	—	

Model (Portable Meters)	Orion Star A300 Series Premium	Orion Star A200 Advanced	Orion Star A100 Series Basic	Orion Intrinsically Safe	
Example Display					
Display Type	Graphic LCD with step-by-step text prompts	Graphic LCD with step-by-step text prompts	Fixed segment LCD with basic text prompts	Fixed segment LCD	
Keypad	Comprehensive with menu-specific function keys	Comprehensive with menu-specific function keys	Simple with highlighted measure key	Simple	
Number of Channels	1 to 3	1	1	100 with date/time (266S, 136S only)	
Datalog	5000 with date/time	1000 with date/time	50	100 with date/time on select meters	
USB Output	1	1	—	—	
RS232 Output	1	1	—	1 (266S, 136S only)	
Methods	10	—	—	—	
Parameter	pH	Star A321 pH meter Star A325 pH/conductivity meter Star A326 pH/RDO/DO meter	Star A221 pH meter	Star A121 pH meter	261S pH meter 266S pH meter
	ISE	Star A324 pH/ISE meter	—	—	—
	Conductivity/TDS/Resistivity/Salinity	Star A322 conductivity meter Star A325 pH/conductivity meter	Star A222 conductivity meter	—	—
	Conductivity/TDS	—	—	Star A122 conductivity meter	131S conductivity meter 136S conductivity meter
	RDO Optical & Polarographic Dissolved Oxygen	Star A323 RDO/DO meter Star A326 pH/RDO/DO meter	Star A223 RDO/DO meter	—	—
	Polarographic Dissolved Oxygen	—	—	Star A123 DO meter	—
	pH/ISE/Conductivity/RDO/DO	Star A329 pH/ISE/conductivity/RDO/DO meter	—	—	—