

Thermo Scientific LP Vortex Mixer

Choose to continuously mix or pulse samples with our low profile vortex mixer.

FEATURES

- Low profile design with small footprint
- Continuous or pulsing functions
- Speed Control: 0 to 300rpm
- Built-in safety features regarding current overload
- Wide variety of optional accessories to accommodate different vessel types





APPLICATIONS

- | Any application that requires vortex mixing in a tube, plate, or glass vessel

11

Ordering Information and Specifications

Cat. No.	88880017	88880018
Plug Configurations		
Speed Range	0-3000rpm	
Operating Mode	Continuous / Touch	
Max Load	0.5 kg / 1.1lbs	
Dimensions (W x D x H) without attachments	154 x 210 x 83mm 6.1 x 8.3 x 3.3in	
Weight	3.1kg / 6.8lbs	
Power	100 – 240V, 50/60Hz	
Certifications	CE	
Warranty	2 Years	

Accessories for the LP Vortex mixer

	Cat. No.
Microplate Tray	88880117
Rubber Platform	88880118
Tube Holder with Insert Retainer	88880119
Tube Holder with Insert Retainer (Φ 10 mm)	88880120
Tube Holder with Insert Retainer (Φ 12 mm)	88880121
Tube Holder with Insert Retainer (Φ 15 mm)	88880122
Tube Holder with Insert Retainer (Φ 20 mm)	88880123
Tube Holder with Insert Retainer (Φ 25 mm)	88880124



Tube Holders