



Thermo Scientific Orion VERSA STAR Advanced Electrochemistry Benchtop Meter & Modules


Don't see the benchtop meter parameter combination you need for your application? No problem – simply order the Orion VERSA STAR meter and choose the modules you need.

With versatility being the name of the game, you can also choose modules to add or change your existing Orion VERSA STAR meter to meet your latest testing needs. Simply swap out the modules when you need to and no need to wait to power down. The meter automatically recognizes the changes and datalogs remain unchanged and stored in the meter.

Each module includes the necessary probe connections for the measurement parameter. The meter includes two Orion Star stirrer probe ports, one RS232 port, two USB ports, power port and grounding port. Included with the meter are the benchtop electrode arm and individual channel covers to protect the channels not in use.

Meter and Module Ordering Information

Cat. No.	Description
VSTAR00	Orion VERSA STAR Benchtop Meter - Electrode arm with redesigned holder - Universal power adapter
VSTAR-PH	Orion VERSA STAR pH/temperature module
VSTAR-ISE	Orion VERSA STAR pH/ISE/temperature module
VSTAR-CND	Orion VERSA STAR conductivity/temperature module
VSTAR-RD	Orion VERSA STAR RDO/DO/temperature module
VSTAR-LR	Orion VERSA STAR pH/temperature with LogR technology module
STARA-BEA	Benchtop electrode arm
STARA-HB	Free-standing, weighted base for use with benchtop electrode arm (sold separately)



No need to turn off your meter to switch modules – meter automatically recognizes module changes.

- pH
- ISE
- ORP
- Cond
- RDO
- DO
- Temp